

Microbiology Tools and Protocols

Measuring tools - Using a Pipette

To set the volume

1. Check the micropipette type.
 - a. For R-1000 and R-2000, turn the volume lock counter-clockwise so the volume setting mechanism is unlocked and free to turn.
 - b. Edvotek does not have a locking mechanism. Skip to step 2.
2. Adjust the volume indicator to the desired quantity by rotating the plunger button. Rotate the plunger button counter-clockwise to increase the volume, and clockwise to decrease the volume.
 - a. For R-1000 and R-2000: The digital volume indicator is read from top to bottom. Red digits = milliliters; Black digits – microliters
 - b. For Edvotek 5-50 μL : The digital volume indicator is read from left to right. Dial reads only μL so set exact volume required.
3. For R-1000 and R-2000 only: Turn the volume lock clockwise to prevent accidental changes to the volume setting.



Aspirating a sample with a Pipette

1. Attach a new tip by inserting the shaft of pipette into the tip.
2. Press the plunger button fully and release your thumb. The piston returns to “home position” and latches automatically, ready to pick up the volume displayed on the volume indicator.
3. Holding the pipette vertically, place the tip into the sample. Press and release the trigger to aspirate the sample. Do not hold the plunger during aspiration.
4. Pause briefly to ensure that the full volume of sample is drawn into the tip.
5. Withdraw the tip from the sample.

Dispensing:

1. Press the plunger slowly to the bottom of the stroke length. Hold for 1-3 seconds.
2. Still holding the plunger, withdraw the tip. Release the plunger. The pipette will return to the home position, ready for the next pipetting cycle.
3. Press the tip ejector button to discard the tip. Use a new tip for each sample to prevent contamination.

Single Channel Digital Pipette

